

Skin prep tissues REF 501

Description

Puradep™ alcohol tissues are designed for cleaning and disinfection of the intact skin prior to intramuscular, intravenous, subcutaneous injections, vaccinations etc.

Packaging information

Packaging	Container
PCS P/C	250 tissues
Selling unit	24 containers
EAN 13 (x01)	8717371750695
EAN 13 (x24)	8717371750015

Specifications

Sterile	No
Material tissue	Cellulose, cotton viscose fibers, polyester fibres and binder
Latex free	Yes



Puradep™ skin prep tissues.

Puradep™ alcohol tissues are designed for cleaning and disinfection of the intact skin prior to intramuscular, intravenous, subcutaneous injections, vaccinations etc. Puradep™ is packed in a handy dispenser container with an easy removable lid. A handy metal wall bracket is available for the Puradep™ container

Instructions for use:

1: Open the lid and remove the alluminium foil. 2: Lead the first tissue from the middle of the roll through the opening of the lid. 3: Close the lid properly.

Puradep™ ingredients (INCI).

0,5% chlorhexidine in 70% ethanol.


Puradep™